

HDF5 MCS Raw Data Definition

Aus Multi Channel Systems Wiki

Inhaltsverzeichnis

- 1 Definition of the HDF5 format for raw data
 - 1.1 Changelog
- 2 Hierarchy
- 3 Root-Folder "/"
- 4 Folder "Data"
- 5 Folder "Recording_x"
- 6 Folder "AnalogStream"
- 7 Sub-folder "Stream_x" of "AnalogStream"
- 8 Folder "FrameStream"
- 9 Subfolder "Stream_x" of "FrameStream"
- 10 Subfolder "FrameDataEntity_x"
- 11 Folder "EventStream"
- 12 Subfolder "Stream_x" of "EventStream"
- 13 Folder "SegmentStream"
- 14 Subfolder "Stream_x" of "SegmentStream"
- 15 **DataSubType-Average**: Subfolder "Stream_x" of "SegmentStream"
- 16 Folder "TimeStampStream"
- 17 Subfolder "Stream_x" of "TimeStampStream"
- 18 Comment

Definition of the HDF5 format for raw data

MCS-HDF5 Protocol Type: **RawData** (Raw-Data protocol)

Protocol Version: **3** based on the definitions of RawDataFileIO in version 10.

All strings are only ASCII-encoded

Changelog

Version 1:

- Initial draft


Version 2:

- New Root-Folder attributes added to detect name and version of the creating application and library

Version 3:

- Data structures for **DataSubType::Average** of **StreamType::Segment** added

Hierarchy



Root-Folder "/"

Contains all information for one experiment - measured data (inside the folder **Data**) and a description (possibly in the future) inside the folder **Experiment/Description/...**

Attributes:

Name	Description	Data Type	MCS-HDF5 Protocol Version
McsHdf5ProtocolType	Type of the used MCS-HDF5 protocol definition (e.g. RawData for the raw data MCS-HDF5 definitions)	[String,Scalar]	1 ≤
McsHdf5ProtocolVersion	Version number of the used MCS-HDF5 protocol	[Integer,Scalar]	1 ≤
GeneratingApplicationName	Name of the application that generated this HDF5 file	[String,Scalar]	2 ≤
GeneratingApplicationVersion	Version of the application that generated this HDF5 file	[String,Scalar]	2 ≤
McsDataToolsVersion	Version of the McsDataTools library that was used by the application to create the HDF5 file	[String,Scalar]	2 ≤

Datasets:

- none

Folder "Data"

Navigation: /Data

Contains all recordings for this experiment.

Attributes:

Name	Description	Data Type
ProgramName	Name of the recording program	[String,Scalar]
ProgramVersion	Version number of the recording program	[String,Scalar]
MeaName	Name of the recorded MEA	[String,Scalar]
MeaLayout	Layout descriptor	[String,Scalar]
MeaSN	Serial number of the MEA	[String,Scalar]
Date	Date of the recording	[String,Scalar]
DateInTicks	Date of the recording in .NET ticks (100 ns)	[Long(64-bit Integer),Scalar]
FileGUID	GUID of the converted raw data file	[String,Scalar]
Comment	Comment	[String,Scalar]

Datasets:

- none

Folder "Recording_x"

Navigation: /Data/Recording_x

Contains all recorded streams for recording x.

Attributes:

Name	Description	Data Type
RecordingID	Recording ID	[Integer(32-bit Integer),Scalar]

RecordingType	Recording type	[String,Scalar]
TimeStamp	Start time of the recording in microseconds	[Long(64-bit Integer),Scalar]
Duration	Total recording duration in microseconds (This duration can differ from the actual duration of the recorded data!!!)	[Long(64-bit Integer),Scalar]
Label	Label	[String,Scalar]
Comment	Comment	[String,Scalar]

Datasets:

- none

Folder "AnalogStream"

Navigation: /Data/Recording_x/AnalogStream

(Organisational) folder for all channel-based streams of this recording

Attributes:

- none

Datasets:

- none

Sub-folder "Stream_x" of "AnalogStream"

Navigation: /Data/Recording_x/AnalogStream/Stream_x

Container for an analog stream

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[Int(32-bit Integer),Scalar]	1 ≤
Label	Label	[String,Scalar]	1 ≤
SourceStreamGUID	GUID of the source streams	[String,Scalar]	1 ≤
StreamGUID	GUID	[String,Scalar]	1 ≤
StreamType	Type of the stream, e.g. Electrode	[String,Scalar]	1 ≤
DataSubType	Sub-type of the analog stream (e.g. Analog)	[String,Scalar]	1 ≤

Datasets:

- Matrix [InfoChannel](#) → $n \times 16$ matrix of describing information vectors for the n channels:
 - **Attributes:** [InfoVersion](#) → Version number of the Info-Objects [Int(32-bit Integer),Scalar]

Name	Description	Data Type	InfoVersion
ChannelID	ID of the channel as given by the recording software	[Int(32-bit Integer),Array(Size 1)]	1 ≤

RowIndex	Row number of this channel inside the ChannelData matrix where the data of this channel is stored	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that this channel belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the channel	[String,Array]	1 ≤
RawDataType	Type of the raw data	[String,Array]	1 ≤
Unit	Physical unit of the measured sensor value	[String,Array]	1 ≤
Exponent	Exponent $n \Rightarrow 10^n$ resp. 10^n in which the channel values magnitude is measured (e.g. k,m, μ ,...)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ADZero	ADC-Step that represents the 0-point of the measuring range of the ADC	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Tick	Sample tick Δ between two sample points of a channel in μs \Rightarrow sampling frequency = $1000000 / \Delta$	[Long(64-bit Integer),Array(Size 1)]	1 ≤
ConversionFactor	Conversion factor for the mapping ADC-Step \Rightarrow measured value	[Long(64-bit Integer),Array(Size 1)]	1 ≤
ADCBits	Number of bits used by the AD-Converter	[Int(32-bit Integer),Array(Size 1)]	1 ≤
HighPassFilterType	Type of the high-pass filter (empty string if not available)	[String,Scalar]	1 ≤
HighPassFilterCutOffFrequency	Cut-off frequency of the high-pass filter ('-1'-String if not available)	[String,Scalar]	1 ≤
HighPassFilterOrder	Order of the high-pass filter (-1 if not available)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
LowPassFilterType	Type of the low-pass filter (empty string if not available)	[String,Scalar]	1 ≤
LowPassFilterCutOffFrequency	Cut-off frequency of the low-pass filter ('-1'-String if not available)	[String,Scalar]	1 ≤
LowPassFilterOrder	Order of the low-pass filter (-1 if not available)	[Int(32-bit Integer),Array(Size 1)]	1 ≤

- 2-dimensional Data-Matrix **ChannelData** \rightarrow Data for sampled channels organized as $n \times m$ matrix \Rightarrow one row per channel and one column per sample time point
 - reconstruct the value of the measured signal:

$$y(\text{channel}, t_{ind}) = (\text{ChannelData}[\text{InfoChannel}[\text{channel}]. \text{RowIndex}, t_{ind}] - \text{ADZero}) * \text{InfoChannel}[\text{channel}]. \text{ConversionFactor} * 10^{\text{InfoChannel}[\text{channel}]. \text{Exponent}}$$
 in $\text{InfoChannel}[\text{channel}]. \text{Unit}$
 - reconstruct the sample time point: $t = t_{ind} * \text{InfoChannel}[\text{channel}]. \text{Tick}$ in μs
- Matrix **ChannelDataTimeStamps** $\rightarrow k \times 3$ matrix of segments where the rows are one segment and the columns are:
 - first column \rightarrow time stamp of the first sample point of the segment
 - second column \rightarrow first index (column) of the segment in **ChannelData**
 - third column \rightarrow last index (column) of the segment in **ChannelData**

Folder "FrameStream"

Navigation: /Data/Recording_x/FrameStream

(Organisational) folder for all frame-based streams of this recording

Attributes:

- none

Datasets:

- none

Subfolder "Stream_x" of "FrameStream"

Navigation: /Data/Recording_x/FrameStream/Stream_x

Folder that contains all Frame-Entities of one Frame-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[Int (32-bit Integer), Scalar]	1 ≤
Label	Label	[String , Scalar]	1 ≤
SourceStreamGUID	GUID of the source stream	[String , Scalar]	1 ≤
StreamGUID	GUID	[String , Scalar]	1 ≤
StreamType	Type of the stream Frame	[String , Scalar]	1 ≤
DataSubType	Sub-type of the event stream (e.g. SpikeTimeStamp)	[String , Scalar]	1 ≤

Datasets:

- Matrix [InfoFrame](#) → n × 24 matrix of describing information vectors for the n Frame-Entities:
 - Attributes:** [InfoVersion](#) → Version number of the Info-Objects [**Int**(32-bit Integer),**Scalar**]

Name	Description	Data Type	InfoVersion
FrameID	ID of the frame entity as given by the recording software	[Int (32-bit Integer), Array (Size 1)]	1 ≤
FrameDataID	ID of the frame entity inside the stream folder that maps this information vector to the entity folder (FrameDataID → subfolder FrameDataEntity_FrameDataID)	[Int (32-bit Integer), Array (Size 1)]	1 ≤
GroupID	ID of the group that this frame entity belongs to	[Int (32-bit Integer), Array (Size 1)]	1 ≤
Label	Label of the entity	[String , Array]	1 ≤
RawDataType	Type of the raw data	[String , Array]	1 ≤
Unit	Physical unit of the measured sensor value	[String , Array]	1 ≤
Exponent	Exponent n ⇒ 1En resp. 10 ⁿ in which the sensor values magnitude is measured (e.g. k,m,μ,...)	[Int (32-bit Integer), Array (Size 1)]	1 ≤
ADZero	ADC-Step that represents the 0-point of the measuring range of the ADC	[Int (32-bit Integer), Array (Size 1)]	1 ≤
ADCBits	Number of bits used by the AD-Converter	[Int (32-bit Integer), Array (Size 1)]	1 ≤
Tick	Sample tick Δ between two frames in μs ⇒ sampling frequency = 1000000 / Δ	[Long (64-bit Integer), Array (Size 1)]	1 ≤
HighPassFilterType	Type of the high-pass filter (empty string if not available)	[String , Scalar]	1 ≤
HighPassFilterCutOffFrequency	Cut-off frequency of the high-pass filter ('-1'-String if not available)	[String , Scalar]	1 ≤
HighPassFilterOrder	Order of the high-pass filter (-1 if not available)	[Int (32-bit Integer), Array (Size 1)]	1 ≤
LowPassFilterType	Type of the low-pass filter (empty string if not available)	[String , Scalar]	1 ≤
LowPassFilterCutOffFrequency	Cut-off frequency of the low-pass filter ('-1'-String if not available)	[String , Scalar]	1 ≤
LowPassFilterOrder	Order of the low-pass filter (-1 if not available)	[Int (32-bit Integer), Array (Size 1)]	1 ≤
SensorSpacing	Distance between adjacent sensors in μm	[Int (32-bit Integer), Array (Size 1)]	1 ≤
FrameLeft	Sensor count of the left edge of the entity frame based on the reference frame	[Int (32-bit Integer), Array (Size 1)]	1 ≤

FrameTop	Sensor count of the top edge of the entity frame based on the reference frame	[Int(32-bit Integer),Array(Size 1)]	1 ≤
FrameRight	Sensor count of the right edge of the entity frame based on the reference frame	[Int(32-bit Integer),Array(Size 1)]	1 ≤
FrameBottom	Sensor count of the bottom edge of the entity frame based on the reference frame	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ReferenceFrameLeft	Sensor count of the left edge of the reference frame (defined by the used sensor array)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ReferenceFrameTop	Sensor count of the left edge of the reference frame (defined by the used sensor array)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ReferenceFrameRight	Sensor count of the left edge of the reference frame (defined by the used sensor array)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
ReferenceFrameBottom	Sensor count of the left edge of the reference frame (defined by the used sensor array)	[Int(32-bit Integer),Array(Size 1)]	1 ≤

Subfolder "FrameDataEntity_x"

Navigation: /Data/Recording_x/FrameStream/Stream_x/FrameDataEntity_x

Contains all datasets of the Frame-Entity x

Datasets:

- Matrix [ConversionFactors](#) → $n \times m$ matrix of conversion factors for the sensor array
- 3-dimensional Data-Cube [FrameData](#) → cube of the frame data organized as one frame to one sample time point ($n \times m$ matrix of sampled signal values per sensor) \times sample time points
 - reconstruct the value of the measured signal: $y = (\text{FrameData}[x,y,t] - \text{ADZero}) * \text{ConversionFactors}[x,y]$
 - reconstruct the sample time point:
- Matrix [FrameDataTimeStamps](#) → $k \times 3$ matrix of segments where the rows are one segment and the columns are:
 - first column → time stamp of the first sample point of the segment
 - second column → first index (z-axis) of the segment in **FrameData**
 - third column → last index (z-axis) of the segment in **FrameData**

Datasets:

- none

Folder "EventStream"

Navigation: /Data/Recording_x/EventStream

(Organisational) folder for all event-based streams of this recording

Attributes:

- none

Datasets:

- none

Subfolder "Stream_x" of "EventStream"

Navigation: /Data/Recording_x/EventStream/Stream_x

Folder that contains all Event-Entities of one Event-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[Int(32-bit Integer),Scalar]	1 ≤
Label	Label	[String,Scalar]	1 ≤
SourceStreamGUID	GUID of the source stream	[String,Scalar]	1 ≤
StreamGUID	GUID of the current stream	[String,Scalar]	1 ≤
StreamType	Type of the stream Event	[String,Scalar]	1 ≤
DataSubType	Sub-type of the event stream (e.g. StgSideband , UserInput , DigitalPort)	[String,Scalar]	1 ≤

Sub-type Description:

- **StgSideband** → The event is associated to a STG sideband change.
- **UserInput** → The event is associated with an user input.
- **DigitalPort** → The event is associated with a digital port change.

Datasets:

- Matrix [InfoEvent](#) → $n \times 7$ matrix of describing information vectors for the n Event-Entities:
 - **Attributes:** [InfoVersion](#) → Version number of the Info-Objects [Int(32-bit Integer),Scalar]

Name	Description	Data Type	InfoVersion
EventID	ID of the event entity	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that the entity belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the entity	[String,Array]	1 ≤
RawDataType	Type of the raw data	[String,Array]	1 ≤
RawDataBytes	Number of bytes of the raw data type	[Int(32-bit Integer),Array(Size 1)]	1 ≤
SourceChannelIDs	Comma separated list of ID's of (source) channel that were involved in the generation of this event	[String,Array]	1 ≤
SourceChannelLabels	Comma separated list of labels of the source channels	[String,Scalar]	1 ≤

- 2-dimensional matrix [EventEntity_x](#) → $5 \times n$ matrix ⇒ n events with describing vector (time stamp of event, duration of event, event info type, info 1, info 2)
 - **Attributes:** Short description of content
 - $t_{\text{event } i} = \text{EventEntity}_x[0, i]$ in μs
 - $\Delta_{\text{event } i} = \text{EventEntity}_x[1, i]$ in μs

Folder "SegmentStream"

Navigation: /Data/Recording_x/SegmentStream

(Organisational) folder for all segment-based streams of this recording. A segment is a cutout of parts of the sampled signal relative to an event, defined by a pre- and post interval.

Attributes:

- none

Datasets:

- none

Subfolder "Stream_x" of "SegmentStream"

Navigation: /Data/Recording_x/SegmentStream/Stream_x

Folder that contains all Segment-Entities of one Segment-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[Int(32-bit Integer),Scalar]	1 ≤
Label	Label	[String,Scalar]	1 ≤
SourceStreamGUID	GUID of the source stream	[String,Scalar]	1 ≤
StreamGUID	GUID of the current stream	[String,Scalar]	1 ≤
StreamType	Type of the stream Segment	[String,Scalar]	1 ≤
DataSubType	Sub-type of the segment stream (e.g. Spike)	[String,Scalar]	1 ≤

Datasets:

- Matrix **InfoSegment** → $n \times 7$ matrix of describing information vectors for the n Segment-Entities:
 - **Attributes:** **InfoVersion** → Version number of the Info-Objects [Int(32-bit Integer),Scalar]

Name	Description	Data Type	InfoVersion
SegmentID	ID of the segment entity	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that the segment entity belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the entity	[String,Array]	1 ≤
PreInterval	Time interval in μs before the segment defining event occurred - definition of the beginning of the segment	[Int(64-bit Integer),Array(Size 1)]	1 ≤
PostInterval	Time interval in μs after the segment defining event occurred - definition of the end of the segment length of the segment = PreInterval + PostInterval in μs	[Int(64-bit Integer),Array(Size 1)]	1 ≤
SegmentType	Type of the segment (e.g. SpikeCutout)	[String,Array]	1 ≤
SourceChannelIDs	Comma separated list of ID's of (source) channels that the segments are taken from → Link to the SourceChannelInfo matrix	[String,Array]	1 ≤

- 2-dimensional matrix **SourceChannelInfo** → $n \times 15$ matrix ⇒ n of describing vectors for the n source channels, the structure is the same as in **ChannelInfo** used in section Sub-folder "Stream_x" of "AnalogStream"
 - **Attributes:** **InfoVersion** → Version number of the Info-Objects [Int(32-bit Integer),Scalar]
- Vector **SegmentData_ts_x** → n time stamps in μs of the event triggering the segment, one for each of the n segments contained by segment entity x
- 2-dimensional matrix or 3-dimensional cube **SegmentData_x** → $k \times n$ matrix (k sample points for one segment, n number of sampled segments) or $k \times m \times n$ cube (k sample points for one segment, m number of segments for one time point/for one multi-segment, n number of sampled multi-segments) of segment data:
 - **Attributes:** **SourceChannelID** → Comma separated list of ID's of (source) channels that the segments are taken from [String,Scalar] (the same as in InfoSegment, repeated for clarification)
 - reconstruct the value of the measured segment signal (only one segment id_{segment} → 2-dimensional matrix $\mathbf{M}[\text{row},\text{col}]$):
 - $t_{\text{ind}}[\text{row}, \text{col}] = \text{SegmentData_ts_x}[\text{col}] + (\text{row} - 1) * \text{tick}_{\text{source-channel}} - \text{PreInterval}$ in μs
 - $y(id_{\text{segment}}, t_{\text{ind}}(\text{row}, \text{col})) = (\text{SegmentData_x}[\text{row}, \text{col}] - \text{ADZero}_{\text{source-channel}}) * \text{ConversionFactor}_{\text{source-channel}} * 10^{\text{Exponent}_{\text{source-channel}}}$ in $\text{InfoChannel}[\text{source-channel}]$. Unit

- reconstruct the value of the measured segment signal (m segments → multi-segments → 3-dimensional cube $\mathbf{M}[\text{row}, \text{col}, z]$):
 - $\text{col} \rightarrow id_{\text{segment}} \rightarrow \text{source-channel}$
 - $t_{\text{ind}}[\text{row}, \text{col}, z] = \text{SegmentData ts } x[z] + (\text{row} - 1) * \text{tick}_{\text{source-channel}[\text{col}]}$ in μs
 - $y(id_{\text{segment}}, t_{\text{ind}}(\text{row}, z)) = (\text{SegmentData } x[\text{row}, \text{col}, z] - \text{ADZero}_{\text{source-channel}[\text{col}]}) * \text{ConversionFactor}_{\text{source-channel}[\text{col}]} * 10^{\text{Exponent}_{\text{source-channel}[\text{col}]}}$ in $\text{InfoChannel}[\text{source-channel}[\text{col}]]$. Unit

DataSubType-Average: Subfolder "Stream_x" of "SegmentStream"

Navigation: /Data/Recording_x/SegmentStream/Stream_x

Folder that contains all Segment-Entities of one Segment-Stream with **DataSubType == Average**:

Attributes: no difference to the standard case above

Datasets:

- Matrix **InfoSegment**: no difference to the standard case above
- Matrix **SourceChannelInfo**: no difference to the standard case above
- $(3 \times n)$ matrix **AverageData_Range_x** → (**start, end, count**) per segment average \times count of segment averages contained by segment entity x. **start** and **end** denote the start and end timestamp in μs of the interval that contains all averaged segments. **count** is the number of averaged segments.
 - **Attributes:** description of the content
- $(2 \times k \times n)$ cube **AverageData_x** → (mean and standard deviation) \times k sample points of the segment \times n number of segment averages
 - **Attributes:** description of the content
 - reconstruct the value of the mean and standard deviation of the average segment (n average segments → 3-dimensional cube $\mathbf{M}[\text{row}, \text{col}, z]$):
 - row: mean → row = 0; StdDev → row = 1
 - col: $t_{\text{ind}}(\text{col}) = (\text{col} - 1) * \text{tick}_{\text{source-channel}} \rightarrow \text{time range } (0, \text{PreInterval}[\text{SegmentID}] + \text{PreInterval}[\text{SegmentID}])$ in μs
 - z: $z = id_{\text{average}}$ (number of average segment)
 - $\text{Mean}(id_{\text{average}}, t_{\text{ind}}(\text{col})) = (\text{AverageData } x[0, \text{col}, id_{\text{average}}] - \text{ADZero}_{\text{source-channel}}) * \text{ConversionFactor}_{\text{source-channel}} * 10^{\text{Exponent}_{\text{source-channel}}}$ in $\text{InfoChannel}_{\text{source-channel}}$. Unit
 - $\text{StdDev}(id_{\text{average}}, t_{\text{ind}}(\text{col})) = \text{AverageData } x[1, \text{col}, id_{\text{average}}] * \text{ConversionFactor}_{\text{source-channel}} * 10^{\text{Exponent}_{\text{source-channel}}}$ in $\text{InfoChannel}_{\text{source-channel}}$. Unit

Folder "TimeStampStream"

Navigation: /Data/Recording_x/TimeStampStream

(Organisational) folder for all TimeStamp-based streams of this recording

Attributes:

- none

Datasets:

- none

Subfolder "Stream_x" of "TimeStampStream"

Navigation: /Data/Recording_x/TimeStampStream/Stream_x

Folder that contains all TimeStamp-Entities of one TimeStamp-Stream:

Attributes:

Name	Description	Data Type	StreamInfoVersion
StreamInfoVersion	Version number of the meta information structure	[Int(32-bit Integer),Scalar]	1 ≤
Label	Label	[String,Scalar]	1 ≤
SourceStreamGUID	GUID of the source stream	[String,Scalar]	1 ≤
StreamGUID	GUID of the current stream	[String,Scalar]	1 ≤
StreamType	Type of the stream TimeStamp	[String,Scalar]	1 ≤
DataSubType	Sub-type of the TimeStamp stream (e.g. NeuralSpike)	[String,Scalar]	1 ≤

Sub-type Description:

- **NeuralSpike** → The entity contains time stamps of neural spikes

Datasets:

- Matrix [InfoTimeStamp](#) → $n \times 7$ matrix of describing information vectors for the n Event-Entities:
 - **Attributes:** [InfoVersion](#) → Version number of the Info-Objects [Int(32-bit Integer),Scalar]

Name	Description	Data Type	InfoVersion
TimeStampEntityID	ID of the event entity	[Int(32-bit Integer),Array(Size 1)]	1 ≤
GroupID	ID of the group that the entity belongs to	[Int(32-bit Integer),Array(Size 1)]	1 ≤
Label	Label of the entity	[String,Array]	1 ≤
Unit	Physical unit of the measured sensor value	[String,Array]	1 ≤
Exponent	Exponent $n \Rightarrow 10^n$ in which the channel values magnitude is measured (e.g. k,m, μ ,...)	[Int(32-bit Integer),Array(Size 1)]	1 ≤
SourceChannelIDs	Comma separated list of ID's of (source) channel that were involved in the generation of this event	[String,Array]	1 ≤
SourceChannelLabels	Comma separated list of labels of the source channels	[String,Scalar]	1 ≤

- Vector [TimeStampEntity_x](#) → n time stamps in μs

Comment

All time-related information except dates (100ns ticks) are given in microsecond ticks!!

-