



MP-RM

Michigan Probe Rod Mount



Example:
Using rod mount
with adapter and
Michigan Probe



It is possible to connect various adapters to the rod mount. Working with different adapters you can add on different probes, for exaxmple MPA8I or MPA32I from Multi Channel Systems GmbH.

Scope of delivery:

- 1 rod mount (stainless steel) for manipulator, length 200mm, diameter 8mm
- 1 grounding cable with multi-disk connector
- 2 screws (M3x6) to fix the adapter for the Michigan Probe