

TC01/2 Temperature controller



Operating temperature	10 °C to 40 °C
Storage temperature	0 °C to 50 °C
Dimensions (W x D x H)	170 mm x 224 mm x 66 mm
Weight	1.5 kg
Main fuse	800 mA slow blow
Supply voltage	85 VAC to 264 VAC @ 47 HZ to 440 Hz
Number of output channels	1 (TC01), 2 (TC02)
Output voltage	max. 15 V
Output current	max. 1 A per channel
Output power	max. 10 W per channel
Resistance	20 Ω
Control range	Ambient temperature (min 10 °C) to 50 °C
Accuracy	0.1 °C
Resolution	0.1 °C