

## ADPT-W2100-HS32-90° for Use with the W2100-System


Adapter to connect a W2100-HS32 Headstage rectangular to the Surface

### Applications

Using two of these rectangular adapters ADPT-W2100-HS32-90° gives the possibility of using two W2100-HS32 headstages at once.

Connect both headstages via adapters vis à vis at each other and record from 64 electrode channels simultaneously.

**ADPT-W2100-HS32-90°**  
**Pin layout of the upright Omnetics Connector**  
Orientate the headstage as shown in the picture.  
Please read also the datasheet W2100-HS32.




ADPT-W2100-HS32-90°

⊗  
E1 to E32  
GND  
REF

Guide post  
Recording electrodes  
Ground  
Reference

ADPT-W2100-HS32-90°  
connected to a W2100-HS32



**ADPT-W2100-HS32-90°**  
**Pin layout for the rectangular Omnetics Connector**  
Please orientate the headstage as shown in the picture.