

Constant Vacuum Pump CVP with Pressure Control

Technical Specification




Constant Vacuum Pump

Top View with Touch Screen



Technical Specification

Operating temperature	10 °C to 40 °C
Storage temperature	10 °C to 50 °C
Maximum pressure	-200 mbar below atmospheric pressure
Resolution	+/- 0.1 mbar
Tension	110 V AC - OR - 230 V AC
Dimensions (H x W x D)	0.47 x 0.15 x 0.23 m
Weight	7.4 kg (without power supply cable)
Sensor type	Differential: Target pressure in relation to actual ambient pressure
Sensor accuracy	+/- 2.5 %
Fuse	1A M
Capacity of the waste bottle	5 l



Rear Panel