

TCW220x100 Warming Plate

Warming Plate for Use with the TC02 Temperature Controller



Warming Plate
TCW220x100

Technical Specification

TCW220x 100	
Dimensions (W x D x H)	220 mm x 100 mm x 10 mm
Weight	600 g
Material	Aluminium, black anodized
Application	Warming plate for in vivo experiments with narcotized rats or mice, and for in vitro experiments, for example, for keeping slices or cell cultures on microelectrode arrays (MEA) warm.
Integrated Sensor	
Integrated sensor:	Type PT 100 1/3 DIN B
Accuracy	± 0.1 °C
Measuring method	Four wire measuring bridge
Resistance	10 -15 W
Heating temperature range	Ambient temperature up to 50 °C
Settling time	± 2 min
Maximal temperature	105 °C
Connection to temperature controller TC02	D-SUD 9 female connector
Cleaning	
Cleaning	Please wipe the warming plate with a soft tissue and detergent or alcohol. Do not immerse the plate into water!